

Prosthetic Feet

Exoskeletal Single Axis Foot

The Sensible Solution

Much like the Endoskeletal Single Axis Foot, the Exoskeletal Single Axis Foot is light-in-weight with an advanced design. This ankle consists of a composite ankle bolted to a wood ankle block, and works well for low activity (U.S. Level 1 or 2) amputees. It features a lifelike cosmesis and interchangeable anterior and posterior bumpers in three durometers.



SPECIFICATIONS

WEIGHT	882 g (1.94 lb)*
PATIENT WEIGHT LIMIT	250 lbs (115 kg)**
HEEL HEIGHT	3/4" (18 mm) 3/8" (10 mm)
AVAILABLE SIZES	22 to 31 cm
TOE RESISTANCE	Low, Regular, High
COSMESIS	Unisex, Lifelike in Buff, Tan or Medium Brown
ACCESSORIES	Bumpers (3 Anterior, 3 Posterior) Exoskeletal Ankles
U.S. WARRANTY	One year and three weeks from purchase
U.S. L-CODE SUGGESTION	5974

- Low Activity (U.S. Levels 1 & 2)
- Lightweight, Advanced Design
- Interchangeable Bumpers in Three Durometers
- Attractive, Lifelike Cosmesis

* Based on size 26 cm foot. Weight includes wood ankle block.

** Body weight plus any loads normally or routinely carried cannot exceed this weight limit.

GENERAL GUIDELINES

The patient's body weight plus any loads normally or routinely carried must not exceed 250 lbs (115 kg). This "modified body weight" is used to determine toe resistance. If the "modified body weight" of the patient exceeds 250 lbs (115 kg), the Pathfinder® II Foot, Magnum Foot or Magnum System is recommended.

Use the patient's "modified body weight" to determine the recommended toe resistance from the chart below. If the patient is of higher activity (U.S. Level 3 or 4) consider a dynamic response foot.

Low (2)	TOE RESISTANCE	
	Regular (3)	High (4)
130 lbs and below	131 - 180 lbs	181 - 250 lbs
59 kg and below	60 - 81 kg	82 - 115 kg



COSMESIS
Medium Brown, Tan or Buff



Prosthetic Feet

ANTERIOR BUMPER SELECTION GUIDELINES

Foot Size	Soft	Medium	Firm
22-23 cm	120 lbs (55 kg) and below	120 - 250 lbs (55-115 kg)	N/A
24-25 cm	N/A	250 lbs (115 kg) and below	N/A
26-27 cm	N/A	185 lbs (84 kg) and below	186-250 lbs (85-115 kg)
28-29 cm	N/A	195 lbs (89 kg) and below	196-250 lbs (90-115 kg)
30-31 cm	N/A	160 lbs (73 kg) and below	161-250 lbs (73-115 kg)

POSTERIOR BUMPER SELECTION GUIDELINES

Foot Size	Soft	Medium	Firm
22 - 23 cm	165 lbs (75 kg) and below	166-250 lbs (76-115 kg)	N/A
24 - 25 cm	175 lbs (79 kg) and below	176-250 lbs (80-115 kg)	N/A
26 - 27 cm	N/A	250 lbs (115 kg) and below	N/A
28 - 29 cm	N/A	195 lbs (89 kg) and below	196-250 lbs (90-115 kg)
30 - 31 cm	N/A	160 lbs (73 kg) and below	161-250 lbs (73-115 kg)

ORDERING GUIDELINES

To order a complete Single Axis Foot Assembly, please select the size and side of the foot, color, toe resistance, heel height, the durometer of the anterior deflection bumper, the durometer of the posterior deflection bumper and that the foot is exoskeletal.

If you wish to order any part of the foot separately, refer to Accessories and Replacement Parts.

	SAA-26LB-2L-MM-EXO	Single Axis Foot SAMPLE Part Number
Name Code	SAA	Single Axis Foot Assembly
Size of Foot	26	22 to 31 cm
Side of Foot	L	L = Left R = Right
Color	B	B = Buff M = Medium Brown T = Tan
Toe Resistance	2	2 = Low 3 = Regular 4 = High
Heel Height	L	L = Low 10 mm (3/8") M = Medium 18 mm (3/4")
Anterior Bumper	M	S = Soft M = Medium F = Firm
Posterior Bumper	M	S = Soft M = Medium F = Firm
Exoskeletal	EXO	Exoskeletal