


Liners & Suspension

Alpha® P-POD™ Liners

Alpha Liners, the fabric-covered gel socket interfaces, have provided a high level of comfort for thousands of adult amputees since 1996. The P-POD™ versions of the Alpha Liner are designed and tested specifically for young amputees and feature the same skin-friendly gel—a thermoplastic elastomer that is infused with mineral oil. The gel cushions the limb and provides protection from shear, impact, and abrasive forces.

Parents who are concerned about their children’s sensitive skin need not worry about the Alpha gel; it has been thoroughly tested at an outside lab to determine its potential to cause a skin reaction. On a scale of 0 to 4, with a score of zero indicating “No reaction” and a score of 4 indicating “Severe reactivity,” the Alpha gel received a score of 0 in every test.

Alpha P-POD Liners are made with our custom-made Alpha MAX fabric. The MAX fabric is three times more abrasion-resistant than many other liners on the market and helps to increase the durability around the trim lines and other common wear areas. This liner is made for those who are highly active and require a heavy-duty liner and have proven to be ideal for children. The Alpha P-POD Locking Liner is designed for use with a shuttle lock and 1/4 x 20 threaded pins.



- **Designed and Specifically Tested for Young Amputees**
- **Made of Skin-friendly, Thermoplastic Elastomeric Gel**
- **Durable MAX Fabric**
- **Comfortable Socket Fit**

SPECIFICATIONS

STYLES	Uniform
COLORS	Buff Medium Brown
THICKNESS	3 mm
SIZES	Extra Small, Extra Small Plus
U.S. CODE SUGGESTION	L-5679 Cushion Liner (bill once per liner) L-5673 Locking Liner (bill once per liner)


For locking liners, Ohio Willow Wood strongly recommends applying Loc-Tite Removable Threadlocker 242 (or equivalent) to the pin threads prior to installing the pin. Failure to ensure that the locking pin is secure in the Locking Liner could result in the pin coming out of the liner, resulting in the potential for injury to the amputee.

Liners & Suspension

ORDERING GUIDELINES: Alpha® P-POD™ Liners

	ALC-5832-E	Alpha P-POD Cushion Liner Sample Part Number
Name Code	ALC	ALC = Alpha Cushion Liner ALL = Alpha Locking Liner
Style	5	5 = Uniform
Fabric	8	8 = MAX Fabric: Buff 9 = MAX Fabric: Medium Brown
Thickness	3	3 = 3 mm
Size	2	2 = Extra Small Plus 9 = Extra Small
Single Liner	E	E = Each

Alpha P-POD Sizing Information



See Alpha Ordering Guidelines for warranty information.