

Appendix

4/2/07

Ohio Willow Wood Fasteners by Type

Part # (if applicable)	Description	Uses
	M4 x 20 mm Socket Head Capscrew	PPF-130000 P-POD 22 mm Tube Clamp Adap. (pylon clamp)
	M5 x 10 mm Socket Head Capscrew	PPF-700250 P-POD Lam. Socket Attach. Block (attach lam. cap)
	M5 x 10 mm Dacromet Socket Head Capscrew	GEO-025 Threaded Clamp Adapter (secure clamp slide) FND-228002 4-Hole Rotatable Rec. (secure clamp slide) FND-228003 4-Hole Rotatable Pyr. (secure clamp slide)
	M6 x 8 mm Socket Head Capscrew	700-AL60X, Lamination Lock (attach Lamination Plate)
	M6 x 10 mm SS Flathead Capscrew	GEO-200, 300, 400, 500, 600 GeoFlex/Lite (attach 4-hole component) FND-135105/06/07 Angle Offset Plate (attach components to Angle plate)
	M6 x 12 mm Flathead Capscrew	GEO-200, 300, 400, 500, 600 GeoFlex/Lite (attach 4-hole component)
	M6 x 10 mm Socket Head Capscrew	700-AIS48X Alpha Lock (attach Forming Plate) 700-AIS10X Locking Lanyard (attach Forming Plate)
	M6 x 12 mm Socket Head Capscrew	700-200/250 Lam. Socket Attachment Block (attach lamination cap)
	M8 x 40 mm Socket Head Capscrew	700-AL60X Lamination Lock (attach Lock Body to Fab. Post)
700-B130	M6 x 12 mm Setscrew	PFDR-XXXX-X Pathfinder II (A/P slide adjustment)
700-B132	M8 x 12 mm Lowhead Capscrew	PFDR-XXXX-X Pathfinder II (toe-out adjustment)
	M8 x 18 mm Socket Head Capscrew	PPF-1540 P-POD Exo Lam. Core (attach lamination cap) PPF-FS006 P-POD Symes Nut
700-B098	M6 x 35 mm SS Flat Head Capscrew	FND-134125/150 Slide Adapter (attach to Alpha Lock or Locking Lanyard)
700-B099	M6 x 35 mm Flat Head Capscrew, Dacromet	MHD-224125/150 Magnum Slide Adapter (attach to Alpha Lock or Locking Lanyard)
	M6 x 40 mm Flat Head Capscrew	FND-995000 Comp. TP Pyramid Adapter (attach to socket)
700-B100	M6 x 30 mm Flat Head Capscrew	700-AIS48X Alpha Lock (attach to socket)
	M6 x 30 mm Flat Head Capscrew, Alloy Steel	700-AIS10X Alpha Lanyard (attach to socket)
700-B102	M6 x 25 mm Flat Head Capscrew	700-375 3/8" 4-Hole Spacer Plate (attach to prosthesis) 700-SP500 Suction Spacer Plate (attach to socket) 700-SP501 Suction Spacer Plate w/Pyr. Rec. (attach to socket) 700-SP502 Suction Spacer Plate w/Pyr (attach to socket) 700-SP503 Suction Spacer Plate w/Slide Adap (attach to socket) 700-SP504 Suction Spacer Plate w/5 mm Offset (attach to socket) 700-SP505 Suction Spacer Plate w/10 mm Offset (attach to socket)
700-B103	M6 x 25 mm Socket Head Capscrew	FND-130018 Alum. 18 mm Adapter (pylon clamp) FND-130030 Alum. 30 mm Adapter (pylon clamp) FND-220030 Titan. 30 mm Adapter (pylon clamp) FND-134130 Alum. 4-Hole 30 mm Adapter (pylon clamp) MHD-002 Magnum 30 mm Adapter (pylon clamp) FND-268000, 001 AK Rotatable 3-Prong Adap., Rec. (clamp) FND-268002 AK Rotatable 3- Prong Base (clamp) STK-222 Stab. Knee Extension Assist (spring adjustment) GEO-100, 200, 500 GeoFlex/GeoLite (pylon clamp)



Appendix

Part # (if applicable)	Description	Uses
700-B104	M6 x 20 mm Flat Head Capscrew	700-SP500 Suction Spacer Plate (attach to socket) 700-SP501 Suction Spacer Plate w/Pyr. Rec. (attach to socket) 700-SP502 Suction Spacer Plate w/Pyr (attach to socket) GEO-025 GeoFlex Threaded Clamp Adapter (attach to socket) FND-228002 4-Hole Rotatable Pyr. Rec. (attach to prosthesis) FND-228003 4-Hole Rotatable Pyr. (attach to prosthesis) PFD-0114 Pathfinder II Through 4-Hole Adap (attach to socket)
700-B105	M6 x 20 mm Flat Head Capscrew, Dacromet	700-300 1/8" 4-Hole Spacer Plate (attach to prosthesis)
700-B106	M6 x 16 mm SS Flat Head Capscrew	700-SP10X Suction Pro (attach 4-hole components) 700-AL60X, Lamination Lock (attach 4-hole components) FND-224120 20mm Offset Pyramid (attach to prosthesis) GEO-200, 300, 400, 500, 600 GeoFlex/Lite (attach 4-hole component) 700-GL481 G-Lock with Pyramid (attach pyramid) FND-135105 Angle Offset Plate (attach Angle Plate to prosthesis) FND-134125/150 Slide Adapter (attach to socket) FND-135105/106/107 Offset Plates (attach plate to prosthesis)
700-B107	M6 x 16 mm Flat Head Capscrew, Dacromet	700-200/250 Lam. Socket Attachment Block (attach 4-hole components) MHD-224125/150 Magnum Slide Adapter (attach to socket)
700-B110	M8 x 25 mm Setscrew	700-AL48X Alpha Lock (attaches Distal Adapter to Fabricating Dome)
700-B114	M10 x 50 mm Socket Head Capscrew	Carbon Copy 2 Men's Foot
700-B116	M10 x 43 mm Socket Head Capscrew	Carbon Copy 2 Women's Foot; Impulse Foot
700-B118	M10 x 38 mm Socket Head Capscrew	SACH Foot
700-B120	3/8" x 2.75" Socket Head Capscrew	English thread for use with all feet
700-B122	M8 x 12 mm Setscrew	Shorter setscrews to use in place of M8-1.25 x 16mm setscrews in cases of extreme alignment
700-B126	M10 x 32 mm Socket Head Capscrew	EW2 Ankle (used with EW2 Foot); Impulse Ankle (used with Impulse
700-B128	M10 x 68 mm Socket Head Capscrew	2LTC-1210 Pyramid Spacer Block (used with Impulse)
700-B134	M6 x 6 mm Setscrew	Trailblazer Foot (attach heel)
700-B200	M8 x 24 mm Socket Head Capscrew	P-POD SACH (with Foot Pyr. Receiver, Symes Nut, or Exo ankle)
700-B201	M8 x 28 mm Socket Head Capscrew	P-POD Impulse (with Foot Pyr. Receiver, Symes Nut, or Exo ankle)
700-B202	M8 x 32 mm Socket Head Capscrew	P-POD SACH (with Foot Pyr. Adapter)
700-B203	M8 x 34 mm Socket Head Capscrew	P-POD Impulse (with Foot Pyr. Adapter)
700-B204	M8 x 55 mm Socket Head Capscrew	PPK-GKIT P-POD Growth Kit
700-B220	M5 x 12 mm Flat Head Capscrew	PPF-227014 P-POD 4-Hole Pyramid (attach to prosthesis) PPF-137004 P-POD 4-Hole Pyramid Rec. (attach to prosthesis) PPF-SP471 P-POD Titanium Suction Pyr. (attach to prosthesis) PPF-700250 P-POD Lam. Socket Attach. Block (attach 4-hole components)
700-B221	M5 x 16 mm Flat Head Capscrew	PPF-994118 P-POD 1/8" Spacer Plate (attach to prosthesis)
700-B222	M5 x 20 mm Flat Head Capscrew	PPF-994114 P-POD 1/4" Spacer Plate (attach to prosthesis)



Appendix

Part # (if applicable)	Description	Uses
700-B223	M6 x 12 mm Setscrew	PPF-137004 P-POD 4-Hole Pyr. Receiver (adjust alignment)
		PPF-130000 P-POD 22 mm Tube Clamp Adap. (adjust alignment)
		PPF-131330 P-POD 300 mm Pylon w/Fixed Adap. (adjust alignment)
		PPF-131340 P-POD 400 mm Pylon w/Fixed Adap. (adjust alignment)
		PPF-138027 P-POD 27 mm Dual Receiver (adjust alignment)
		PPF-138051 P-POD 51 mm Dual Receiver (adjust alignment)
		PPF-1209 P-POD Alum. Foot Pyr. Rec. (adjust alignment)
700-FS002	M8 x 16 mm Setscrew	FND-130018 Alum. 18 mm Adapter (adjust alignment)
		FND-130030 Alum. 30 mm Adapter (adjust alignment)
		FND-220030 Titan. 30 mm Adapter (adjust alignment)
		FND-131330 Al. Fixed Adap. w/30 mm tube (adjust alignment)
		FND-221330 Ti. Fixed Adap. w/30 mm tube (adjust alignment)
		FND-139930 Al. Fixed Adap. w/23 mm rod (adjust alignment)
		FND-229930 Ti. Fixed Adap. w/23 mm rod (adjust alignment)
		FND-228002 4-Hole Rotatable Pyr. Rec. (adjust alignment)
		FND-268000 AK Rotatable 3-Prong Pyr. Rec. (adjust alignment)
		FND-268004 Threaded Rotatable Pyr. Rec. Insert (adjust alignment)
		FND-228075 75 mm Dual Receiver (adjust alignment)
		FND-228060 60 mm Dual Receiver (adjust alignment)
		FND-228045 45 mm Dual Receiver (adjust alignment)
		FND-228032 32 mm Dual Receiver (adjust alignment)
		MHD-002 Magnum 30 mm Adapter (adjust alignment)
		MHD-PR-4 Magnum Titan. 4-Hole Pyr. Receiver (adjust alignment)
PFD-0116 Pathfinder II Pyramid Receiver (adjust alignment)		

