

Joint Solutions

Stabilized Knee

The Stabilized Knee is the perfect choice for less active transfemoral amputees. This weight-activated stance control knee features a stainless steel wave spring extension assist — allowing for a smoother gait, while aiding in noise reduction. Its strong, lightweight titanium design incorporates this weight-activated stance control feature, and both distal and proximal connectors.



| Product No. | Description |
|-------------|--|
| STK-221 | Titanium Stabilized Knee (Modular) |
| STK-222 | Titanium Stabilized Knee (Modular), with extension assist |
| STK-223 | Stainless Steel Stabilized Knee (Modular) |
| STK-224 | Stainless Steel Stabilized Knee (Modular), with extension assist |

- U.S. Activity Levels 1 & 2
- Available in Titanium or Stainless Steel
- Extension Assist Design Featuring Stainless Steel Spring Extension Aid
- Weight Activated Stance Control Feature
- Proximal and Distal Pyramid Connections

SPECIFICATIONS

WEIGHT

| | |
|-----------------|--|
| Titanium | Knee 357 g + Extension Assist 90 g = 447 g |
| Stainless Steel | Knee 525 g + Extension Assist 90 g = 619 g |


PATIENT WEIGHT LIMIT 250 lbs/115 kg*

U.S. L-CODE SUGGESTION Knee: 5812; Extension Assist: 5850

U.S. WARRANTY If a knee is returned to Ohio Willow Wood for repair within the one-year warranty period, a loaner will be supplied during the time of repair service. If a knee is returned outside the warranty period, a loaner will be available for a nominal fee and service quotes will be given prior to repair. During the warranty period, the \$50 Loaner fee will be waived.

* Body weight plus any loads normally and routinely carried cannot exceed this weight limit.

Stabilized Knee Accessories



Replacement Knee Covers, Spring Adjustment Clamps, and Extension Assists are available.

| Product No. | Description |
|-------------|---------------------------|
| 1 STK-003 | Spring Adjustment Clamp |
| 2 STK-004 | Protective Cover |
| 3 STK-100 | Extension Assist Assembly |